
Cadwell® Sleepmate

Instructions for Use (IFU)

Thermocouple

Contents	Page
English	4
Dansk.....	6
Deutsch.....	8
Español	10
Suomi	12
Français.....	14
Italiano.....	16
Nederlands.....	18
Norsk	20
Português.....	22
Svenska.....	24

Intended use

The Cadwell Sleepmate™ Reusable Airflow Sensor is to detect patient nasal and/or oral breathing during sleep research and clinical studies of sleep.

Instruction for use

The sensor should be installed so that the probes are exposed to the maximum available airflow. The sensor probes are flexible to allow each probe to be bent **gently** to insure that they are in the path of airflow. The probes should not touch the skin as contact may attenuate the signal amplitude. All cannula style airflow sensors are installed onto the patient in a manner similar to that of an oxygen cannula. **The sensors are delicate and must be handled carefully.** Proper handling and cleaning will increase durability and longevity in the life of the sensor. When bending the prongs to fit the patient, never bend at a sharp angle, but in a smooth round curve. The prong should be bent a short distance away from the center bar. The assembly should be carefully secured using both the sliding cinch tubing and surgical tape.

Connect the sensor to the appropriate channel input of your recorder. No external power source is required.

Ensure the sensor is securely attached to the patient as described and you have clear, reliable signals prior to lights out.

Cleaning

This sensor was designed to be reusable and can be cleaned by rubbing with any cleansing agent that is not corrosive toward plastic. **The entire assembly must be completely dry before using again to insure patient safety.** For sterilization implement a standard procedure for gas sterilization. Do not autoclave this product. Use the molded ends of the connectors when removing the connector from the junction box.

Technical Specifications

Description	Properties
Lead Wire	Two conductors white lead wire, approx. 200 cm (80") long.
Sensor	Flexible, rubber coated probes for oral (one probe) and nasal (two probes) airflow.
Connectors	Two 1.5 mm (0.1") safety pins. Some recorders may use Lemo or other connectors.
Filter Settings	Low Hz = 0.1 Hz; High Hz = 15 Hz; Gain: Output range is +/- 100 MV. Sensitivity is dependent on variables such as probe placement, patient airflow, ambient room temperature, CPAP pressure, etc.

Caution Statement:

US Federal Law restricts this device to sale by or on the order of a physician. This product is for diagnostic purposes only and is not intended to be used as an apnea monitor or in a life supporting or sustaining capacity.

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©2021 Cadwell Industries, Inc. All rights reserved. PN 302993-639-5 Rev. 00
Cadwell Industries Inc.
1 (800) 245-3001 U.S.A. toll-free | +1 (509) 735-6481
This document contains trademarks that belong to Cadwell Industries, Inc.