
Cadwell® Sleepmate

Instructions for Use (IFU)

Snore Sensor

Contents**Page**

English	4
Dansk	6
Deutsch	8
Español	10
Suomi	12
Français	14
Italiano	16
Nederlands	18
Norsk	20
Português	22
Svenska	24

Intended use

The Cadwell Sleepmate™ Snore Sensor or Microphone is intended to detect patient snoring during sleep research and clinical studies of sleep.

Instruction for use

When positioning the snoring microphone, place it on a location of the neck that will be least offensive to the patient. Do not place the microphone or sensor directly on the larynx. When positioning the snoring sensor it is best to have the patient simulate a snore while you place your fingers gently on the neck. Try to select a spot that presents the most “vibration” associated with the sound. **It is important that the sensor adhere closely to the skin to insure that the vibrations of each snore are transmitted to the sensor.** Place the snore sensor so that the side with the raised bubble is next to the skin.

Attach using tape intended for facial use. Using an additional piece of tape about 15 cm (six inches) away from the microphone or sensor can be useful for strain relief.

Connect the sensor to the appropriate channel input of your recorder. No battery or interface is needed.

Ensure the sensor is securely attached to the patient as described and you have clear, reliable signals prior to lights out.

Cleaning

You can clean the device by rubbing with any cleansing agent that is not corrosive toward plastic.

The entire assembly must be completely dry before using again to insure patient safety. For sterilization implement a standard procedure for gas sterilization. Do not autoclave this product. Use the molded ends of the jackpins when removing the connector from the junction box.

Technical Specifications

The snoring microphone and snoring sensor will work on virtually all standard physiographs or computer recording equipment.

Description	Properties
Lead wire	Two conductor white, approx. 200 cm (80").
Connectors	1.5 mm (0.1") safety pins or special connectors such as a Lemo on some models.
Microphone:	
Sensitivity level	2 - 75 μ V / mm
Filter settings	High Pass 70 Hz or 120 Hz, Low pass 5 Hz or 0.03 sec
Microphone dimensions	Round disc diameter = 16 mm (0.6")
Sensor:	
Filter Settings	High Pass: 50 Hz; Low Pass: 250 Hz.
Signal Output	Approx. 1 mV for loud snoring
Output Impedance	50 k Ω
Sensor Dimensions	Round flexible silicon rubber disc (diameter = 20 mm; thickness = 5 mm)

Caution Statement:

US Federal Law restricts this device to sale by or on the order of a physician. This product is for diagnostic purposes only and is not intended to be used as an apnea monitor or in a life supporting or sustaining capacity

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©2021 Cadwell Industries, Inc. All rights reserved. PN 302994-639 Rev. 00
Cadwell Industries Inc.
1 (800) 245-3001 U.S.A. toll-free | +1 (509) 735-6481
This document contains trademarks that belong to Cadwell Industries, Inc.