
Cadwell® Sleepmate

Instructions for Use (IFU)

Piezo Respiratory Effort Sensor

Contents**Page**

English	4
Dansk	6
Deutsch	8
Español	11
Suomi	14
Français	16
Italiano	19
Nederlands	22
Norsk	25
Português	27
Svenska	30

Intended use

The Cadwell Sleepmate™ Respiratory Effort Sensor is intended to detect respiratory effort from chest or abdomen movement during sleep research and clinical studies of sleep.

Instruction for use

The system takes advantage of the piezo electric sensor's ability to produce a small voltage when stress is applied. On each breath the sensor is stressed by the expansion of the thorax or abdomen and the voltage produced is amplified by the physiograph to produce a clean, reliable respiratory effort signal.

This respiratory system consists of two parts: (1) a sensor which is encased in molded box or glued between elastic strips with a black plastic strip extending on both sides, and (2) a latex free belting system for attaching the sensor to the patient. Please note that we offer three styles of sensors. This User Guide is written for all three styles.

The Belts are intended to be worn over patient nightclothes.

Attach the band to the buckle or loop on one end of the sensor. Place the sensor on the patient. Wrap the band around the patient and insert the remaining free end of the band into the buckle or loop on the other side of the sensor. Adjust for a comfortable fit. You should be able to insert 2 fingers between the band and the patient. Be sure to grasp the buckle or loop of the sensor unit in one hand during adjustment to minimize the tension that is applied directly to the sensor. NO interface or external power source is necessary because the sensors produce their own signal.

When placing the thoracic band, it is frequently best NOT to place the sensor in the center of the breasts. Moving the sensor slightly off center can avoid artifact associated with obese or large breasted patients. Running the lead wire under the band provides additional strain relief.

Connect the sensor to the appropriate channel input of your recorder.

Ensure the sensor is securely attached to the patient as described and you have clear, reliable signals prior to lights out.

Cleaning

You can clean the system by rubbing with any cleansing agent that is not corrosive toward plastic. The belting system can be machine washed in gentle cycle or hand washed. **The entire assembly must be completely dry before using again to insure patient safety.** For sterilization implement a standard procedure for gas sterilization. Do not autoclave this product. Use the molded ends of the jackpins when removing the belt from the junction box.

Technical Specifications

Description	Properties
Lead wire	Two conductor white, approx. 200 cm (80").
Connectors	Two 1.5 mm Safety pins, some models may use special connectors.
Signal output	Output signal is typically 1 to 3 mV @ 200 g, 0.5 Hz.
Belting system	The sensors use velstretch bands of varying lengths depending on patient needs. Disposable bands are also available.
Filter/Time constant	0.1 Hz / One second or longer. The use of a shorter time constant may attenuate waveforms.
High Frequency Filter	15 Hz
Sensitivity	Approximately 30 μ V/mm. Adjustment of the sensitivity may be required. Response is dependent on a variety of variables such as sensor application or patient effort.

Caution Statement:

US Federal Law restricts this device to sale by or on the order of a physician. This product is for diagnostic purposes only and is not intended to be used as an apnea monitor or in a life supporting or sustaining capacity

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©2021 Cadwell Industries, Inc. All rights reserved. PN 302995-639-2 Rev. 00
Cadwell Industries Inc.
1 (800) 245-3001 U.S.A. toll-free | +1 (509) 735-6481
This document contains trademarks that belong to Cadwell Industries, Inc.