




Cadwell® Sleepmate

Instructions for Use (IFU)

RIPmate™ Respiratory Effort

Contents**Page**

English	4
Dansk	9
Deutsch	13
Español	18
Suomi	23
Français	28
Italiano	33
Nederlands	38
Norsk	43
Português	48
Svenska	53

Intended Use

The Cadwell Sleepmate™ respiratory effort sensor is intended to detect respiratory effort from chest or abdomen movement during sleep research and clinical studies of sleep. This Belt is intended for use in supervised medical settings.

Instruction for use

a How to Adjust the Size

The shortest length of the belt is when the sliders are positioned at each end of the belt, as illustrated in figure: *Shortest position*

For ease of adjustment, it is best to keep the sliders centered (equal distance from the end of belts to each slider) in the belt to ensure ease of adjustment.

b Lengthening the Belt

1. Move the left slider towards the center of the belt. A loop will be formed as shown in the picture.
2. Holding the left end of the belt, pull the right end to tighten the loop.

c Shortening the Belt

To shorten the belt, move the sliders apart from each other towards the belt's ends.

d Patient Attachment

The Belts are intended to be worn over patient nightclothes.

The belt should fit snugly around patient without being uncomfortable.

1. Adjust belt as described above for proper fit.
2. Wrap 1 belt around patient Abdomen and snap buckles together.
3. Plug the 1.5mm touch proofed plug at the end of the wire set into the belt connectors near each buckle; plug the other end, with keyhole connector into the RIPmate™ Abdomen module.
4. Plug the module output into the appropriate bipolar input of your polygraph. The red touch proof is plugged into (+) input, and the black into (-) input.
5. Repeat steps 1-4 for the Thorax belt, wrapping this belt around the patient chest.

Note: Do not use two abdomen or two thorax modules in the same recording. This will cause interference between the two sensors, resulting in bad signals.

Warnings

- Do not use damaged belts or cables.
- Do not use the device in an MRI environment. Do not use the device in an explosive environment, i.e. in the presence of flammable liquids or gases, such as anesthetic mixture with air or with oxygen or nitrous oxide
- The RIPmate™ systems are only to be connected to an electrically isolated input. Do not plug cables into electrical outlets as this could result in serious electrical shock.
- The device contains a battery and must therefore be disposed of properly. Local, state or national laws may prohibit disposal of batteries in ordinary trash. Contact your local waste authority for information regarding available recycling and disposal options.
- The belts shall be worn over night clothes.
- The RIPmate™ system is non-repairable by the user and the user's service personnel. Contact Sleepmate in case of device failure.
- Portable and mobile RF communications can affect the performance of the device.
- No user serviceable parts inside. Serviced by Sleepmate and authorized parties only. Warranty void if opened.

Cleaning

No part of the RIPmate™ system requires sterilization.

Sensor

The sensor may be wiped clean with hospital grade cleaner that **is not corrosive to plastic or metal**, and then dried with a clean, dry cloth. Do not immerse the sensor in liquid and avoid contact of the cleaning solution with the connectors.

Belts

The belts may be hand washed in warm (not hot) solution of household laundry detergent, and then air dried. Care should be taken to fully dry safety jacks after cleaning.

Cables

The cables may be wiped clean with a moist cloth using hospital grade cleaner.

Do not immerse the cables in liquid and avoid contact of the cleaning solution with the connectors.

Storage

Proper storage of the RIPmate™ system will increase life of the product. To protect the belt and the cable from damage, follow these guidelines for storage:

- Do not fold the belt tightly together (see below) as this could damage the wire weaved in the belt. Hanging on a hook is the recommended storage method.
- Do not wrap the cable around the sensor as it may cause the cable to break.
- **To save battery power, the module MUST be disconnected from the belt when not in use. When disconnecting cable from belt, grasp both belt connector & cable connector securely and pull out. Do not pull on wire.**

Technical Specifications

Description	Properties
Physical Properties	
Belt length	Adjustable 24"-79" (60-200cm)
Belt material	Elastic band with Insulated wire
Cable length	80" (205cm)
Cable material	PVC jacketed; ZIP-style cord
Module	Box: ABSCable: PVC jacketed; Zip Style cord
Battery	3V Lithium (integrated in the device)
Battery Life	The battery life is estimated to be 3000 hours or 375 eight hour studies. The battery is activated once the sensor module is connected to the belt.
Environmental Specifications	
Temperature	Operation : 40°F to 120°F (+5°C to +50°C) Storage: 0°F to 120°F (-20°C to +50°C)
Relative Humidity	Operation : 15-95% (non-condensing) Storage: 10-95% (non-condensing)
Pressure	Withstands atmospheric pressures from 7.3 psi to 29 psi
Output Specifications	
Output Signal	Maximum signal amplitude: ±5mV
Frequency Range	0.2 to 3 Hz
Sensitivity	Approximately 50µV/mm
Polygraph Settings	
Sensitivity	Adjustment of the sensitivity up or down is typically required. Response is dependent upon variables such as sensor application and patient effort.
Low Frequency Filter / Time Constant	0.16 Hz (or 1 second or longer) Shorter time constants or higher low frequency filter settings will significantly attenuate waveforms
High Frequency Filter	35 Hz.

**Caution Statement:**

US Federal Law restricts this device to sale by or on the order of a physician. This product is for diagnostic purposes only and is not intended to be used as an apnea monitor or in a life supporting or sustaining capacity.

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©2021 Cadwell Industries, Inc. All rights reserved. PN 302995-639-4 Rev. 00
Cadwell Industries Inc.
1 (800) 245-3001 U.S.A. toll-free | +1 (509) 735-6481
This document contains trademarks that belong to Cadwell Industries, Inc.