

	Safety Data Sheet / Scheda di sicurezza	
	Product/Prodotto : SAC2	
Date/Data 27/03/2019	Rev. 0	

SAFETY DATA SHEET SAC2

1. Product and Firm identification

1.1 *Preparation Identification* SAC2

1.2 *Preparation Used* SAC2 is an adhesive EEG electrode cream to use with surface electrodes. SAC2 is intended to be applied on intact skin.

1.3 *Firm Information* Spes Medica Srl
Via Buccari 21, 16153 Genova (GE) – Italy
e-mail: info@spesmedica.com
web: www.spesmedica.com

1.4 *Emergency Phone number* +39 010 39 03 43

2. Hazards Identification

2.1 This product is not considered dangerous either for human or environment.

3. Composition/Ingredients information

Components: Water, Talc, Celite, Glycerol, CarboxyMethylCellulose, Sodium chloride, Potassium chloride, Phenoxyethanol, Ethylhexylglycerin

4. Medical Aid Measures

4.1 *Skin* Generally the product does not irritate the skin. The product is intended for use on intact skin.

4.2 *Eyes* Wash immediately and continuously with flowing water for 15 minutes. Otherwise, contact the doctor.

4.3 *Ingestion* Contact the doctor.

4.4 *Inhalation* N.A.